




just smile

Dental Implant



A minimally invasive solution in the atrophic posterior maxilla

JD PTERYGO


A MINIMALLY INVASIVE APPROACH

The JDPterygo Ø 4mm implant is manufactured in titanium grade 4, with a 3.5mm anodised implant collar.

It is designed to obtain an optimal primary stability, achieving a perfect anchorage in the pterygoid bone.

A minimally invasive approach is possible also with computer guided surgery.





NEW Ø 3.3mm SOLUTION

FOR EXTREME ATROPHY

LESS THAN 6mm OF BONE

The new JD Pterygo Ø 3.3mm implant is manufactured in titanium grade 4, with a 3.5mm anodised implant collar.

It is especially designed for extreme atrophy less than 6mm of bone between the distal wall of the sinus and the end of the tuberosity.

With a simplified drilling protocol with only two drilling steps, you will easily obtain a perfect anchorage in the pterygoid bone.



JD PTERYGO


PRODUCT SPECIFICATIONS

The JDP terygo dental implants are available in the diameters and lengths shown in the following table:





JDPterygo Ø 3.3mm connection is compatible with JDEvolution S



JDPterygo Ø 4mm connection is compatible with JDEvolution Plus+

IMPLANT DIAMETER	TIP DIAMETER	LENGTH			
Ø 3,3	1,75	13	15	18	20
Ø 4,0	2,4	13	15	18	20

Note: All measurements in mm

PRODUCT CATALOGUE

Implants:

a	2	ာ
U	J	,0

 PT40130S:
 JDPterygo Ø 3.3 L 13

 PT40150S:
 JDPterygo Ø 3.3 L 15

 PT40180S:
 JDPterygo Ø 3.3 L 18

 PT4020OS:
 JDPterygo Ø 3.3 L 20

Ø 4,0

PT40130: JDPterygo Ø 4.0 L 13
PT40150: JDPterygo Ø 4.0 L 15
PT40180: JDPterygo Ø 4.0 L 18
PT40200: JDPterygo Ø 4.0 L 20



Drills for non guided protocol:

JDDRPT20Drill Ø 2.0 JDPterygoJDDRPT24Drill Ø 2.4 JDPterygoJDDR101Drill Ø 2.8 JDPterygoJDDRPT32Drill Ø 3.2 JDPterygo



Guided drills for JDPterygo guided protocol:

JDGD20-130	Guided Drill Ø 2.0 L 13.0
JDGD24-130	Guided Drill Ø 2.4 L 13.0
JDGD20-150	Guided Drill Ø 2.0 L 15.0
JDGD24-150	Guided Drill Ø 2.4 L 15.0
JDGD20-180	Guided Drill Ø 2.0 L 18.0
JDGD24-180	Guided Drill Ø 2.4 L 18.0
JDGD20-200	Guided Drill Ø 2.0 L 20.0
JDGD24-200	Guided Drill Ø 2.4 L 20.0




Implant and prosthetic drivers:

JDID100 JD Implant Driver - JDPterygo

EVSDPF60: Prosthetic Driver for JDTorque L 60 JDEvolution Plus

JDID101 JD Implant Driver - JDPterygo One EVSDPF60 Prosthetic Driver for JDTorque L 60



JD Pterygo Surgical Kit:

JDKIT01JD Surgical Kit - JDPterygo Drills KitJDKIT02JD Guided Surgery Kit Extra DrillsJDKIT04JD Surgical Kit - JDPterygo Kit

Prosthetic Solutions:

JDPterygo Ø 3.3mm connection is compatible with JDEvolution S, so please refer to JDEvolution S catalogue to choose the most suitable component. JDPterygo Ø 4mm connection is compatible with JDEvolution Plus, so please refer to JDEvolution Plus catalogue to choose the most suitable component.

JD PTERYGO

SITE PREPARATION DRILLING SEQUENCE

JDP terygo implant site preparation can be performed following both non-guided and computer guided surgical protocol. It is recommended to adhere to the indications of the following drilling sequence to ensure optimal primary stability of the implant.

JD PTERYGO Ø 3.3

· Non-guided protocol



1. Start the osteotomy using JDPterygo drill Ø 2.0mm at the same implant length to be inserted



2. Complete the osteotomy with JDPterygo drill Ø 2.8mm at the entrance for 6mm



The JDPterygo Drills Kit consists of four longer drills with two laser marks, ideal for working in the posterior region of the maxilla



The JD Pterygo Kit consists of four longer drills, two screwdrivers and two implant drivers.

Computer guided protocol



- 1. Start the osteotomy using the Guided Drill Ø 2.0mm at the same implant length.
- 2. Complete the osteotomy with Guided Drill Ø 2.8mm for 6mm.

JD PTERYGO Ø 4.0

Non-guided protocol



1. Start the osteotomy using JDPterygo drill Ø 2.0mm at the same implant length to be inserted.



2. Continue using JDPterygo drill Ø 2.4mm at the same implant length to be inserted.



3. Complete the osteotomy with JDPterygo drill Ø 3.2mm at the entrance for 6mm.

Computer guided protocol



- . Start the osteotomy using the first Guided Drill JDGD24-060 Ø 2.4 x L 6mm*
- 2. Continue using Guided Drill Ø 2.4 at the same implant length to be inserted (Choose among Guided drills L 13-15-18-20mm)
- 3 Continue with the Guided Drill Ø 2 8mm for 8mm*
- 4. Complete the osteotomy with Guided Drill Ø 3.2mm for 6mm*

*the marked drills are not present in this brochure.
Please refer to JD Guided Surgery Kit brochure available at www.jdentalcare.com



www.jdentalcare.com